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# Instructions for use Non-Sterile Orthopedic Instruments

Valid from: 2017 – 01 - 01

Version: 1.0

TREU Instrumente GmbH  
take-off Gewerbepark 130-132  
D-78579 Neuhausen ob Eck  
Phone: +49 7467 94 76 73 0  
Fax: +49 7467 94 76 73 - 50

## Products

These instructions for use are valid for non-sterile orthopedic instruments.

## Important note

*Read these instructions for use carefully before every application and store them in a location which is easily accessible to users or professional technicians.*

*Read warning notices designated with this symbol carefully. Inappropriate use of these products can severely injure patients, users, or third parties.*

## 1 Area of Application

The products listed above may only be used by appropriately qualified and trained personnel. Products may only be used in a sterile environment.

## 2 Precautionary Measures and Warnings

- Do not use steel wool or abrasive cleaning agents.
- Do not use any cleaning solutions with high iodine or chlorine concentrations.
- Do not place contaminated or used Treu medical products into a case for cleaning in a cleaning / disinfecting machine. Contaminated Treu products must be prepared separately from screens and cases. Cases are designed as organizing containers for steam sterilization, as storage containers for medical products, and as organizing containers during an operation.
- Clean long, narrow cannulations, blind holes, and complicated components with special care.
- Clean all products thoroughly.
- Treu instruments must be sterilized in their final packaging before use.
- Sterilization parameters apply exclusively to components that have received adequate cleaning in advance.
- The parameters listed apply exclusively to properly installed, maintained, and calibrated treatment systems that fulfill the requirements of standards ISO 15883 and ISO 17665.
- Patients at risk due to Creutzfeldt-Jakob disease (CJD) and associated infections

must be operated on with disposable instruments. Instruments used to operate on a patient suspected of having CJD or shown to have the disease must be disposed of after the operation, and / or the currently valid national recommendations must be observed.

– For further information, see applicable national laws and guidelines. Internal guidelines and process instructions within the clinic, and the recommendations and instructions of the cleaning and disinfecting agent and clinical processing system manufacturers must also be observed.

## 3 Limits of Clinical Processing

– Repeated processing cycles involving ultrasound cleaning, mechanical cleaning, and steam sterilization have only a minimal effect on Treu surgical instruments.

– In general, signs of wear and use should determine when a product has reached the end of its service life. Signs of damage and wear include, but are not limited to: Corrosion (rust, pitting), discoloration, deep scratches, flaking, erosion, and cracks. Instruments that are not working properly, or defective instruments subject to excessive use, and instruments with illegible markings or missing or removed (abraded) part numbers must not be used.

## 4 Handling Instruments in the Operating Room

Frequently clean instruments of blood and/or tissue debris during operations to prevent them from drying on.

- Rinse cannulated instruments thoroughly with sterile or processed water to prevent contamination and / or tissue debris from drying onto the internal surfaces.
- To protect personnel and the environment, separate contaminated instruments from non-contaminated items.
- Cover contaminated instruments with a cloth dipped in sterile or processed water to prevent blood and / or tissue debris from drying onto them.

## 5 Condition on Delivery

Medical products are delivered non-sterile, and must be processed and sterilized by the user in accordance with the following instructions before the first use, and then before every subsequent use.

## 6 Preparation

### Additional Information

- If municipal water supplies are used, they must correspond to GUIDELINE 98/83/EC OF THE COUNCIL of November 3rd, 1998 on the Quality of Water for Human Consumption.
- The cleaning and disinfecting agents used for the validation are listed in

these processing instructions. If the processor uses an alternative cleaning agent or disinfecting agent (RKI or VAH listed), they are responsible for the consequences of this decision.

- Re-assemble disassembled products before sterilization.

### Transport

Products must be disposed of without drying out immediately after use. This means that products must be transported wet in a closed container from their location of use to processing, so as to ensure nothing dries onto the products. Transport contaminated items separate from non-contaminated products to protect against contamination.

### Preparation for Decontamination

If possible, products must be disassembled before completing the following processing steps or opened during further processing steps. Ensure no areas are left unwashed. Products must be processed in suitable screen baskets or drainage containers (select size according to the specific product). Products must be attached inside the cleaning basket with a minimum spacing between them. Do not allow products to overlap in order to avoid damaging them during the cleaning process.

### Manual Pre-Cleaning

#### A) Brushing

If possible, disassemble instruments and clean with a brush under cold water <40°C until all visible contamination is removed.

#### B) Cleaning with Water Jet

Rinse instruments with Lumen or Luer Lock fittings with a water jet (4 bar) > 10 sec. Use the water jet to rinse out cavities, joints, gaps, and holes.

#### C) Ultrasound

If pre-cleaning with a brush and rinsing with a water jet still leaves behind visible contamination, complete a pre-cleaning with ultrasound technology.

Dip instruments in an ultrasound bath filled with water <40°C, 0.5 % alkaline detergent (Neodisher FA, Dr. Weigert) for 5 minutes and clean for 10 minutes. Rinse instruments with a water jet (4 bar) > 10 sec.

### Cleaning / Disinfection

#### Automatic cleaning / disinfection process:

(Washing machine, Washer - Disinfector G 7735 CD (Miele):

Step 1: Pre-clean for 1 minute(s) with cold municipal water, drinking water quality <40°C

Step 2: Drain water

Step 3: Pre-clean for 3 minute(s) with cold municipal water, drinking water quality <40°C

Step 4: Drain water

Step 5: Clean for 5 minute(s) at 55°C±5°C with 0.5 % alkaline cleaning agent (Neodisher FA, Dr. Weigert).

Step 6: Drain water

Step 7: Neutralize for 3 minute(s) with cold municipal water, drinking water quality <40°C

# Instructions for use Non-Sterile Orthopedic Instruments

Valid from: 2017 – 01 - 01

Version: 1.0

## Step 8: Drain water

**Step 9:** Rinse for 2 minute(s) with cold municipal water, drinking water quality <40°C

Observe specialized instructions from the cleaning equipment manufacturer.

### Automatic Disinfection

Automatic thermal disinfection in cleaning and disinfecting device, observing national requirements, at A<sub>0</sub>- value 3000:

>5 minutes at 92°C±2°C  
with demineralized water.

### Automatic Drying

Automatic drying in accordance with the automatic drying process of the cleaning and disinfecting machine for at least 30 minutes (at 60°C±5°C in the washing chamber). Complete the following manual drying as needed with a lint-free cloth, and blow out lumen using sterile- oil-free compressed air.

### Sterilization

Sterilize products using the fractional pre-vacuum process (in accordance with DIN EN ISO 17665-1), observing applicable national requirements. Products must be sterilized in suitable sterilization packaging in accordance with ISO 11607 and EN 868.

### EU Standard

Sterilization must be completed using a fractional pre-vacuum process with the following parameters:

134°C

≥5 minutes wait time

3 pre-vacuum cycles

Drying in vacuum for at least 20 minutes

### US Standard

Sterilization must be completed using a fractional pre-vacuum process with the following parameters:

132°C

≥4 minutes wait time

3 pre-vacuum cycles

Drying in vacuum for at least 20 minutes

Observe the autoclave manufacturer's instructions for use and recommended guidelines for the maximum sterilized material load. The autoclave must be installed, maintained, validated, and calibrated in accordance with requirements.

### Additional Information

The processor is responsible for ensuring that the preparation completed, with the equipment, materials, and personnel employed in the preparation facility, achieve the desired results. In general, validation and routine monitoring of the process and the equipment used will be required.

## 7 Inspection

Inspect Treu instruments after preparation and before sterilization for the following aspects:

- Cleanliness
- Damages, including but not limited to signs of corrosion (rust, pitting), discoloration, deep scratches, flaking, cracks, and erosion
- Proper function, including but not limited to the sharpness of cutting tools, bending of flexible instruments, motion of hinges/joints/case locks and movable components such as handles, ratchets, and couplers
- Missing or removed (abraded) part numbers
- Instruments that are not working properly, or defective instruments subject to excessive use, and instruments with illegible markings or missing or removed (abraded) part numbers must not be used.

Inspect instruments to ensure surfaces are in good condition, and to ensure they are functional and assembled correctly. Do not use severely damaged instruments, or instruments with illegible markings, signs of corrosion, or blunt cutting edges.

Lubricate instruments with movable components such as hinges, joint connections, spring-mounted ball bearings, and threaded parts. We recommend using only the recommended specialized oil to maintain and care for TREU instruments (e.g. STERILIT oil spray JG 600 or JG 598 care oil, physiologically harmless in accordance with DAB (German Pharmacopoeia) and sec. 31 LMBG (German Act on Foodstuffs and Goods in Daily Use).

### Service and Repair

Do not complete any repairs or modifications to the product yourself. Only authorized manufacturer personnel are responsible for such repairs. If you have any complaints, claims, or suggestions related to our products, please contact us.

### Return Transport

Defective or non-conforming products must have gone through full reprocessing before they are returned for repair / service.

## 8 Storage and Disposal

Packaged sterile products should be stored in a dry, clean environment, protected from direct sunlight, pests, extreme temperatures, and humidity.

Products, packaging materials, and accessories must be disposed of in accordance with nationally valid provisions and laws. The manufacturer provides no specific instructions on disposal.

## 9 Liability and Warranty

TREU Instrumente GmbH as manufacturer is not liable for subsequent damages that occur due to improper use or handling. This also applies to repairs or modifications to the product completed by personnel not authorized by the manufacturer. These exclusions of liability also apply to warranted services.

## 10 Description of Symbols

|   |                              |
|---|------------------------------|
|  | Warning!                     |
|  | Observe instructions for use |
|  | Item number                  |
|  | Batch designation            |
|  | CE mark                      |
|  | Non-sterile                  |
|  | Manufacturer                 |

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