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# Instructions for use for External Finger Fixators and Mandible Fixator

CE  
0483

## **Intended Use / Indication**

External finger and mandible fixators are designed for the following indications on phalanges and metacarpal bones in the hand, and to fixate simple mandible fractures:

- Closed comminuted fractures
- Open fractures
- Dislocated joint fractures which can be repositioned through ligamentotaxis
- Bone, joint, and soft tissue infections
- Complex soft tissue injuries
- Bone defects caused by trauma or tumor resection

### **Note**

Suitable fixators may only be selected and used by medical professionals with surgical training, since the fixator must be selected based on the specific bone defect, degree of activity, and any comorbidities. Care must be provided in accordance with common traumatology standards, and the implant selected must be appropriate for the type of fracture.

- The doctor is obligated to explain the advantages and disadvantages of the fixator to the patient.

## **Contraindication**

- 1.) Physical activity associated with major physical concussions that could expose the implant to shocks and/or excessive loads.
- 2.) Allergies to a material component.
- 3.) Patients with psychological or neurological problems who are not capable of following the instructions for post-operative treatment.
- 4.) Severe osteoporosis (Singh index 1-3), persistent diabetes mellitus, and HIV infection may represent a counterindication.
- 5.) Manifest abuse of medications, drugs, or alcohol
- 6.) Insufficient or poor quality hard or soft tissues.

## **Possible Side Effects**

- 1.) Nerve and/or vascular lesions
- 2.) Deep or surface infections in the pin tract, osteomyelitis, and chronic secretion at the pin penetration points after screws are removed.
- 3.) Joint contractures, subluxation, luxation, and/or restricted movement.
- 4.) Delayed fracture healing or development of pseudarthrosis.
- 5.) Re-fracture after fixator removal.
- 6.) Loose and/or broken implants
- 7.) Persisting or chronic recurring complications requiring treatment.
- 8.) Second intervention to exchange a module or the full assembly.
- 9.) Allergy to implant steel.
- 10.) Tissue necrosis.
- 11.) Hematoma and poor wound healing.
- 12.) Risks related to anesthesia
- 13.) Chronic pain

## **Notes and Precautionary Measures**

- 1.) If the fracture is fresh, we do not recommend compression.
- 2.) The length of the implant, its thread and diameter must be selected depending on the specific dimensions of the bone.
- 3.) All products must be inspected to ensure they are ready for use before being applied. Defective or damaged parts may not be used.
- 10.) The fixator must have a minimum spacing to the patient's skin, in consideration of post-operative swelling. A certain spacing is also necessary to clean and care for the areas where screws enter the skin.
- 4.) Inspect the fixator assembly regularly to ensure it is stable, and tighten fixator screws if necessary.
- 5.) Careful care of the screw penetration points will avoid pin complications.
- 6.) Patients should receive exact instructions for how to care for their fixator pins.
- 7.) Patients must be instructed that they can turn to the treating physician at any time with questions, complications, or if something is unclear.
- 8.) The patient's healing process must be monitored clinically and radiologically on a regular basis, and fixators must be adjusted / corrected as needed.
- 9.) For patients having distraction completed to elongate a bone, monitor the distraction speed and make adjustments depending on how well the bone is regenerating. The typical distraction speed is 1mm per day.
- 10.) The fixator must be removed after healing is ensured through clinical and radiological means.

## **Material**

Treu Finger and Mandible fixators consist of various components and modules made of implant steel, the titanium alloy Ti6AL4V, and stainless steel.

# Instructions for use for External Finger Fixators and Mandible Fixator

CE  
0483

## Compatibility

For metallurgic, mechanical, and construction-related reasons, implants by different manufacturers and from different materials may never be combined. Material information is provided in the product catalog or on the product labels. Before beginning treatment, ensure that required instruments are available and that they are suitable for combination with our implants.

**Treu can accept no liability if its products are used in combination with third party products.**

The following implants must be used:

Kirschner wires:  
Diameters  $\varnothing$  1.6 and 1.8 mm;  
Length in mm: 70

## Magnetism

Although these fixators are made of stainless and non-magnetic steel or titanium alloys, the magnetic fields used in MRI examinations could accidentally cause them to move or heat up.

Therefore, using magnetic fields is not recommended.

## Notes and Warnings

**Fixators are not designed for re-use. If they are placed under additional strain after they have already been used once, we cannot ensure that they will stand up to this new load.**

**Re-using fixators can lead to severe injuries or even death of the patient!**

A label with a batch number is attached to the fixator packaging. In order to ensure full traceability of the fixator, the physician must enter this number on the patient's operation report.

## General Principles of Hygiene and Preparation

Treu fixators are delivered exclusively in a NON STERILE state. The label includes a symbol to this effect.

### Information for fixators delivered in a non-sterile state:

- Brand new fixators labeled as non-sterile must be prepared before initial use. Transport packaging, protective caps, etc. are not suitable for sterilization.
- Only use permitted materials (RKI, DGHN, VHA)
- Both alkali and PH neutral cleaning agents can be used
- Water quality must correspond to DIN EN 285 attachment B
- Only validated cleaning / disinfection / sterilization processes with sufficient device and product specificity must be applied
- Manufacturer information and recommendations must be observed
- Due to the product design and materials used, we are not able to establish any defined limit of max. processing cycles. This medical product must be treated carefully, and its service life must be determined based on its function.
- Defective products must have gone through full reprocessing before they are returned for repair.

## Preparation for Cleaning

Fixators must be separated and cleaned and sterilized individually, if possible, before the following preparatory steps

### Cleaning

#### Manual cleaning:

- Products must be rinsed under running municipal water ( $>40^{\circ}\text{C}$ ) until all visible contamination is removed
- If necessary, a soft brush should be used to remove the visible contamination
- Immerse instruments into an enzymatic cleaning agent (if an ultrasound bath is used, ultrasound processes of 3 minutes and an ultrasound frequency of 35 KHz are effective).
- Follow the instructions of the cleaning agent manufacturer
- Rinse the instrument under running municipal water ( $>40^{\circ}\text{C}$ )

#### Machine cleaning

place the products in a sterilization tray on the rack trolley and start the cleaning process

- Rinse 1 minute(s) with cold water
- Empty
- Rinse 3 minute(s) with cold water
- Empty
- Wash for 5 minutes at  $55^{\circ}\text{C}$  with 0.5% alkaline cleaner
- Empty
- 3 minutes neutralization with warm tap water ( $>40^{\circ}\text{C}$ ) and neutralizer
- Empty
- 2 minutes intermediate flush with warm tap water ( $>40^{\circ}\text{C}$ )
- Empty

Specialized requirements from the cleaning machine manufacturer must be observed

# Instructions for use for External Finger Fixators and Mandible Fixator

CE  
0483

## **Disinfection**

### **Manual Disinfection**

Dip products into cold water for at least 5 min. Clean the products under cold water with a soft brush until no more residue is visible. Pressure wash threads with a water jet for at least 10 sec. (pulse). Place products into an ultrasound bath at 40°C with 0.5 % enzymatic cleaner for 15 minutes and treat with ultrasound.

After chemical disinfection and cleaning, products must be thoroughly rinsed with clear, flowing water. Any remaining residual contamination must be removed manually (do not use a metal brush or abrasive agent!) To avoid water spots, we recommend a final rinse with completely demineralized water. Then the fixators must be dried immediately.

Follow the manufacturer instructions for the disinfecting cleaner. Ensure that the cleaning agent can reach all interior and exterior parts of the fixator. Clean the product after the contact time with distilled water to remove the cleaning agent.

### **Machine Disinfection**

Complete machine thermal disinfection in consideration of national requirements with regard to A0 value (see ISO 15883).

## **Drying**

Dry the outside of the products by putting them through the dryer cycle of the cleaning / disinfection device. If necessary, you may also clean the products manually with a lint-free cloth. Dry cavities in products with sterile compressed air.

## **Inspection, Maintenance, and Testing**

Fixators must be inspected for discoloration, notches, cracks, and other damage that could be caused by improper sterilization and/or storage before use.

If one of the features above is found on the fixator, it may not be used or inserted before being inspected by us.

Visual inspection for cleanliness.

## **Packaging**

Standard instrument packaging for sterilization in accordance with DIN EN ISO 11607 and DIN EN 868 et seqq.

# Instructions for use for External Finger Fixators and Mandible Fixator

CE  
0483

## Sterilization

- Steam sterilization is recommended
- Other sterilization processes and flash sterilization are not permitted.
- Use steam sterilization at 134° C, with a time of 5 minutes and a pressure of 2.3 bar (according to DIN EN ISO 17665-1; ANSI AAMI ISO 11134)

Cleaning and sterilization must only be completed by qualified professional technicians. Furthermore, improper preparation could cause the fixator to be damaged, making it unsuitable for use.

See validation report MDT no. 04m068 dated July 28th, 2004.

## Storage

Sterilized instruments must be stored in a dry, clean, and dust-free environment at moderate temperatures.

## Disposal

After completing disinfection, defective or explanted fixators must be disposed of professionally.

## Graphic Symbols / Labels

Symbols used for labeling in accordance with DIN EN 980 or DIN EN ISO 15223-1 have the following meanings:

Not for re-use

Manufacturer symbol

Observe instructions for use

Observe accompanying documents!

CE mark with identification number of indicated location

Batch number

Item number

Indicates non-sterile product

Quantity

Only for the (treating) physician

Treu Instrumente GmbH  
Gewerbepark take-off 130-132  
78597 Neuhausen ob Eck / Germany